

Environmental Toxicology Laboratory
Toxicology Centre
University of Saskatchewan

STANDARD OPERATING PROCEDURE

UofS-ETL-EDNA-22

Metabarcoding of Fish Communities

Version 1, April 2018

Yuwei Xie, Ph.D. and John P. Giesy, Ph.D., FRSC, FSETAC

Supported through:
Toxicology Centre and
Department of Veterinary Biomedical Sciences

Correspondence to:
Environmental Toxicology Laboratory
Toxicology Centre
44 Campus Drive,
Saskatoon, Saskatchewan, S7N 5B3
Canada

Phone: (306) 966-5062; 966-2096
Facsimile: (306) 966-4796

APPROVAL PAGE

Revisions to an existing SOP, addition of a SOP change form, or preparation of a new SOP must be reviewed, approved, and signed by the following:

Authored By: Yuwei Xie and John P. Giesy Date: 04/12/2018

Supervisor Review By: John P. Giesy Date: 04/30/2018

Reviewed By: _____ Date: _____
(QA Coordinator)

DEFINITIONS AND ACRONYMS

ETL	Environmental Toxicology Laboratory (University of Saskatchewan)
DQO	Data Quality Objective
DHSE	Department of Health Safety and Environment (University of Saskatchewan)
QA	Quality Assurance
QAPP	Quality Assurance Project Plan
SOP	Standard Operating Procedure
GWF	Global Water Futures
eDNA	Environmental DNA
eRNA	Environmental RNA
PCR	Polymerase chain reaction

TABLE OF CONTENTS

Section	Heading	Page
1.0	PURPOSE	5
2.0	SCOPE AND APPLICATION	5
3.0	SAFETY CONSIDERATIONS	5
	3.1 PERSONAL PROTECTIVE EQUIPMENT	5
4.0	EQUIPMENT, MATERIALS, AND REAGENTS	5
5.0	METHOD, PROCEDURES, AND REQUIREMENTS	6
	5.1 Mobilization and Training	6
	5.2 Decontamination of instruments and laboratory areas before and after isolation	7
	5.3 Prepare primer stock solution and work solution	7
	5.4 PCR amplification for metabarcoding of multiple communities	7
6.0	RECORDS, DOCUMENTATION, AND QC REQUIREMENTS	9
7.0	RESPONSIBILITIES	9
8.0	REFERENCES	10
	APPENDIX I PRIMER SET	11

1.0 PURPOSE

This SOP is developed for environmental DNA (eDNA) metabarcoding of fish communities.

2.0 SCOPE AND APPLICATION

This SOP applies to the ETL for eDNA metabarcoding from the Global Water Futures (GWF) program titled “Next generation solutions to ensure healthy water resources for future generations” (eDNA project).

3.0 SAFETY CONSIDERATIONS

Biosafety training, laboratory training and medical monitoring requirements are described in the Health and Safety Plan for eDNA project of the GWF program. In addition, there are various safety concerns regarding the preparing, storage, and disposal of sample containers at the ETL. If biological spill happens, the appropriate personnel will be immediately notified and the Department of Health, Safety and Environment (DHSE) will be called in order to assess the hazard. DHSE will also be contacted in the case of chemical spills and will be responsible for the disposal of hazardous wastes.

3.1 PERSONAL PROTECTIVE EQUIPMENT

Personal protective equipment (PPE), consisting of nitrile gloves will be worn at all times when handling samples.

4.0 EQUIPMENT, MATERIALS, AND REAGENTS

- Freezer (-20 °C, -80 °C)
- Permanent markers
- Kim wipes
- Nitrile gloves
- Absolute ethanol
- Centrifuge for 0.2 mL 96-well PCR plate ($\leq 4000 \times g$)
- Microcentrifuge (20,000 $\times g$)
- PCR thermal cycler

- 15 mL sterile Graduated Conical-Bottom Tubes with rack
- 1.5 mL sterile microcentrifuge tubes
- Adjustable Pipettes (1000 µl, 200 µl and 10 µl)
- Disposable aerosol filter pipette tips (1000 µl, 200 µl and 10 µl)
- PCR PRE-station
- Qubit V4
- Qubit® dsDNA HS Assay kit
- Qubit® assay tubes (500 tubes, Life Technologies, Cat. no. Q32856)
- Spray bottle with Invitrogen™ RNase Away™ Decontamination Reagent
- Spray bottle with 70% ethanol
- Spray bottle with 10% bleach
- Spray bottle with di-Water
- Thermo Scientific Phusion Green High-Fidelity DNA Polymerase

5.0 METHOD, PROCEDURES, AND REQUIREMENTS

The following procedures describe methods for metabarcoding of fish communities in a stepwise fashion and explain the reasoning behind the techniques. Detailed below are the preparatory procedures, molecular techniques, and procedures for sample documentation and preservation.

5.1 Mobilization and Training

This section describes requirements for mobilization for the necessary lab-work, as well as the necessary safety training. The objective of this section is to ensure that all of the necessary preparatory work has been conducted to enable the successful completion of the overall project. Mobilization for the necessary lab-work entails maintenance and decontamination of equipment and biosafety training in accordance with the site Health and Safety Plan (HASP).

Note: The following precautions must be observed to minimize the risk of cross contamination when setting up a PCR reaction.

- PCR should ideally be performed in a dedicated ‘clean’ area which is free from other work involving DNA.
- All stocks of pipette tips, microcentrifuge tubes, etc. to be used should be DNase, RNase and nuclease-free, and stored in a dust free environment. The use of filter tips is advisable.
- Pipettes designated solely for PCR should be used to set up the reaction if available.
- Wear gloves at all times.
- Prepare the PCR master mix in a “Master Mix laboratory” which is DNA/RNA isolation and gel gel-electrophoresis free lab.

- Template DNA/cDNA should be added in the "DNA extraction laboratory".

5.2 Decontamination of instruments and laboratory areas before and after isolation

- A. Wipe down surfaces of all points of contact, including bench tops, pipettes, and fridge/freezer handles, using 10% bleach. Let bleach remain on surfaces for 15 minutes.
- B. After 15 minutes, use a DI water-dampened paper towel to remove bleach residue.
- C. Use a 70% alcohol dampened paper towel to help quickly dry the surfaces. **Note: For PCR PRE-station, DO NOT use buffers containing ethanol to clean the clear acrylic surface.**
- D. If necessary, repeat the cleaning procedure from step B.
- E. For biosafety cabinet and PCR PRE-station, turn on the UV-light for at least 15 min.
- F. Decontaminate laboratory equipment during work using RNase AWAY™ Surface Decontaminant (Thermo Scientific™)
- G. Wipe down surfaces of bench tops, pipettes, and racks using RNase AWAY™ Surface Decontaminant, then wipe dry with dust-free paper or rinse clean.

5.3 Prepare primer stock solution and work solution

- A. Preparing primer stock solution for long-term storage
Reconstitute dual-tagged forward and reverse primers to 100 µM using 1 x TE Buffer (DNase and DNA free). Carefully transfer the primer stock solutions into the 96-well plate (primer stock solution plate).
- B. Preparing primer working solution for PCR setup
Reconstitute dual-tagged primers to 10 µM using DNase and RNase free sterile deionized water (diH₂O). Add 90 µL diH₂O into each well of a new 200 µL 96-well plate (primer working solution plate), then transfer 10 µL primer stock solution into the primer working solution (10 µM) plate. Seal the working solution plate with a film, mix the solution and centrifuge at 1500 rpm for 2 min. Transfer 20 µL working solution from each well of both the forward and reverse primer working solution plate into a new 96-well plate (primer mixture working solution plate, 40 µL, totally 10 µM for each dual-tagged primer pair, 5 µM for each primer) to prepare a primer mixture working solution. Dual-tagged primer pairs are listed in Appendix I at the end of this document.

5.4 PCR amplification for metabarcoding of multiple communities

Take care to avoid carry over between tubes or contamination of stock.

PCR amplification is conducted in triplicate to reduce the bias introduced by the PCR reaction.

Note: Include the following control reactions for each DNA template.

- Positive DNA control if appropriate
- Negative DNA control if appropriate
- No template control using sterile water instead of DNA

Metabarcoding of multiple barcode genes is utilized to profile the composition of fish communities, including mitochondrial 12s rRNA gene (Primer set: teleoF and teleoR; length: ~100 bp), mitochondrial cytB gene (Primer set: L14841F and 15149R; length: ~470 bp) and mitochondrial COI gene (Primer set: mlCOIintF and jgHCO2198R; length: ~330 bp). The primer sets are dual-tagged to distinguish individual samples.

A. Set up master mix of PCR reaction

a. For mitochondrial 12s rRNA gene and mitochondrial cytB gene

<i>Component</i>	<i>20-μL rxn</i>
diH ₂ O	11.8 μ L
5X Phusion Green HF Buffer	4.0 μ L
10 mM dNTP mix	0.4 μ L
DMSO	0.6 μ L
10 μ M primer mix	2.0 μ L
Phusion Hot Start II DNA Polymerase	0.2 μ L
Template DNA or controls	1.0 μ L
Total	20 μ L

b. For mitochondrial COI gene

<i>Component</i>	<i>20-μL rxn</i>
diH ₂ O	15.0 μ L
10X High Fidelity PCR Buffer	4.0 μ L
50 mM MgSO ₄	0.8 μ L
10 mM dNTP mix	0.4 μ L
10 μ M primer mix	0.8 μ L
Platinum™ II Taq HS HF DNA Polymerase	0.08 μ L
Template DNA or controls	1.0 μ L
Total	20 μ L

B. PCR amplification

Sit the PCR tubes or 96-well plate on a thermocycler, run the amplification program below.

a. For mitochondrial 12s rRNA gene and mitochondrial cytB gene

	<i>Step</i>	<i>Temperature</i>	<i>Time</i>
A	Initial Denaturation	98 °C	30 s
B	Denaturation	98 °C	30 s
C	Annealing	58 °C	30 s

	<i>Step</i>	<i>Temperature</i>	<i>Time</i>
D	Extension	72 °C	30 s
E	Cycling	From B to D	40 cycles
F	Final extension	72 °C	10 min
G	Hold	4 °C	∞

b. For mitochondrial COI gene

	<i>Step</i>	<i>Temperature</i>	<i>Time</i>
A	Initial Denaturation	94 °C	2 min
B	Denaturation	94 °C	15 s
C	Annealing	65 °C	30 s
D	Extension	68 °C	30 s
E	Touchdown	From B to D, -1 °C per cycle	16 cycles
F	Denaturation	94 °C	15 s
G	Annealing	49 °C	30 s
H	Extension	68 °C	30 s
I	Cycling	From F to G	25 cycles
J	Final extension	72 °C	10 min
K	Hold	4 °C	∞

- C. Quality check the PCR products using agarose gel-electrophoresis;
- D. Check and clean the PCR products or store it at -20 °C.

6.0 RECORDS, DOCUMENTATION, AND QC REQUIREMENTS

Preserve PCR products under -20 °C. Sample management will also follow a scripted procedure (ETL-SOP-EDNA-09: Management of eDNA sample: Receiving, Preservation, Storage, Documentation, Decontamination, and Disposal) to ensure samples will not be lost or destroyed. The results of quality control of eDNA samples should be archived properly (hard-copy and electric copy).

7.0 RESPONSIBILITIES

Project Director — Will oversee and approve all project activities.

Project Manager — Will oversee and approve all project activities; review QA reports; approve final project QA needs; authorize necessary actions and adjustments to accomplish program QA objectives; and act as liaison between agencies, field staff.

Quality Assurance (QA) Manager — Will oversee all QA activities to ensure compliance with contract specifications; initiate audits on work completed by project personnel and subcontractors, including analytical laboratories and independent data validation contractors; review program QA activities, quality problems, and quality-related requests. In response to field and analytical findings, this person will approve the corrective actions. This person will report quality non-conformances to the Project Manager and review all pertinent portions of the deliverables before they are transmitted to ensure conformance with QA/QC procedures and quality work product.

Data Manager — Will oversee data management for this project. This person is responsible for the structure, organization, format, implementation, and operation of the project database.

Field Team Leader — Will oversee field activities and supervise the field crews. This person will ensure that proper sample collection, preservation, storage, transport, and COC QC procedures are followed. This person will inform the Project QA Manager when field problems occur, and will communicate and document corrective actions taken. The Field Team Leader will discuss field activities with the Project Manager.

Laboratory Project Manager — The laboratory project manager for this project is the person responsible for assuring that the analysis of all samples is performed in accordance with the QAPP and the laboratory's quality assurance manual. In addition, the Laboratory Project Manager performs the final laboratory review of project data packages for completeness and compliance with project requirements.

8.0 REFERENCES

- Invitrogen ezDNase Product Information Sheet. Pub. No. MAN0015899. Rev. A.0.
Invitrogen SuperScript™ IV First-Strand Synthesis System User Guide. Pub. no. MAN0013442 Rev. B.0.
Thermo Scientific Phusion Green High-Fidelity DNA Polymerase User Guide. Pub. No. MAN0012395. Rev. Date 23 September 2016 (Rev. B.00).
Thermo Scientific. Qubit™ 4 Fluorometer manual. Pub. No. MAN0017209. Rev. B.0.
Thermo Scientific. Qubit® dsDNA HS Assay Kits manual. Pub. No. MAN0002326. Rev. B.0.

APPENDIX I PRIMER SET

L14841F and 15149R PLAT-I

	1	2	3	4	5	6	7	8	9	10	11	12
A	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A
B	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B
C	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C
D	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D
E	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E
F	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F
G	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G
H	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H

L14841F and 15149R PLAT-II

	1	2	3	4	5	6	7	8	9	10	11	12
A	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B
B	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C
C	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D
D	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E
E	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F
F	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G
G	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H
H	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A

teleoF and teleoR PLAT-I

	1	2	3	4	5	6	7	8	9	10	11	12
A	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A
B	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B
C	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C
D	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D
E	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E
F	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F
G	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G
H	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H

teleoF and teleoR PLAT-II

	1	2	3	4	5	6	7	8	9	10	11	12
A	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B
B	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C
C	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D
D	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E
E	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F
F	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G
G	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H
H	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A

mIColintF and jgHCO2198R PLAT-I

	1	2	3	4	5	6	7	8	9	10	11	12
A	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A
B	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B
C	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C
D	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D
E	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E
F	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F
G	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G
H	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H

mIColintF and jgHCO2198R PLAT-II

	1	2	3	4	5	6	7	8	9	10	11	12
A	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B	R.B
B	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C	R.C
C	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D	R.D
D	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E	R.E
E	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F	R.F
F	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G	R.G
G	F01	F02	F03	F04	F05	F06	F07	F08	F09	F10	F11	F12
	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H	R.H
H	F13	F14	F15	F16	F17	F18	F19	F20	F21	F22	F23	F24
	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A	R.A